

BOSON

In Vitro Diagnostic Devices

CATALOG

(Rapid Test)

Introduction

Boson Biotech Co.,Ltd., as a specialist of in vitro diagnostic kits field, we develop and manufacture high-quality point of care and other immunoassay kits for world-wide market.

Our 5,000 square meter facility is operated strictly under ISO 13485:2008 and GMP guidelines. We combine the technological, material and labor resources of China and US to develop, manufacture and supply high quality in vitro diagnostic products efficiently and cost effectively. Currently, we are offering product lines of rapid tests, CLIA, ELISA and raw materials. These product lines provide the immunoassays in various formats to detect cardiac markers, drugs of abuse, fertility hormones, infectious diseases, tumor markers and animal diseases. Many of our products have been approved by Chinese SFDA.

Boson Biotech's mission is to provide the affordable high quality products to help fighting diseases and illicit substance abuse.



Professional Manufacturer of Rapid Diagnostic Tests

Boson Biotech Co., Ltd. is a leading manufacture of Medical Diagnostic Test in China. The company was established in 2001, located in China's special economic zone, Xiamen. With highly experienced staff, ISO13485 compliance, Boson offers premium quality products or OEM, directly from the production facility. The company is efficiently structured to provide lasting healthcare solutions to domestic and international market.



OEM Capabilities for Customers Supply and Manufacturing

Boson also provides a number of customized services to professional distributors and partnering affiliates. The production fully supports branding and product customization and can deliver finished products within two weeks after ordering.

Strongly Integrated Quality Pricing & Service Policies

Every aspect of the company's operations is strictly measured according to Quality, Pricing & Service policies.

Boson provides an exciting opportunity for distribution of quality products at every competitive price because of:

- Integrated production facility, quality control and R&D team, International marketing.
- Most of all, its highly reputable staff with expertise in all aspects of its operation.




HIV 1/2

One Step HIV 1/2 Test (Whole Blood / Serum / Plasma)

Rapid HIV I&II Test is a single use, rapid test device, double antigen sandwich immunoassay for qualitative detection of antibodies to Human Immunodeficiency Viruses or HIV I&II in serum, plasma or whole blood . It is intended for medical institution as the aid for the diagnosis of possible HIV infection. The product may be used for the quick screening of blood donors and blood products.

Sensitivity: 99.6%, Specificity: 99.7%.

This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1N01S2	HIV I&II Test Strip	50 or 100
1N01C2	HIV I&II Test Card	25 or 40
1N01C2-2	HIV I&II Triline Test Card	25 or 40
1N01C2-3	HIV Self Test Kit	1

HCV

One Step HCV Antibody Test (Whole Blood / Serum / Plasma)

Rapid HCV Antibody Test is a single use, rapid test device, double antigen sandwich immunoassay for qualitative detection of antibodies against hepatitis C virus (HCV Ab) in human serum, plasma or whole blood. It is intended for medical institution as an aid for diagnosis and management of patients related to infection with hepatitis C as well for primary screening of blood from volunteer donors on the spot.

Sensitivity: 97.6%, Specificity: 99.5%.

This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1H08S2	HCV Antibody Test Strip	50 or 100
1H08C2	HCV Antibody Test Card	25 or 40
1H08C2-2	HCV Antibody Self Test Kit	1

HBV (HBsAg, HBsAb, HBeAg, HBeAb, HBcAb)

One Step HBV Test (Whole Blood / Serum / Plasma)

Rapid HBV Test is a rapid immunochromatography assay for the qualitative detection of the markers of Hepatitis B virus including Hepatitis B virus surface antigen (HBsAg), Hepatitis B virus surface antibody (anti-HBs/HBsAb), Hepatitis B virus envelope antigen (HBeAg), Hepatitis B virus envelope antibody (HBeAb/anti-HBe) and Hepatitis B virus core antibody (anti-HBc/HBcAb) in human serum, plasma or whole blood specimens. It is intended for use in medical institution as an aid for diagnosis and management of patients related to infection with hepatitis B as well for screening of blood donors or blood products.

HBsAg: Sensitivity: 99.6%, Specificity: 99.6%.

HBsAb: Sensitivity: 97.6%, Specificity: 98.2%.

HBeAg: Sensitivity: 99.4%, Specificity: 99.6%.

HBeAb: Sensitivity: 96.9%, Specificity: 98.8%.

HBcAb: Sensitivity: 97.1%, Specificity: 97.8%.

This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1H01S2	HBsAg Test Strip	50 or 100
1H01C2	HBsAg Test Card	25 or 40
1H02S2	HBsAb Test Strip	50 or 100
1H02C2	HBsAb Test Card	25 or 40
1H03S2	HBeAg Test Strip	50 or 100
1H03C2	HBeAg Test Card	25 or 50
1H04S2	HBeAb Test Strip	50 or 100
1H04C2	HBeAb Test Card	25 or 40
1H05S2	HBcAb Test Strip	50 or 100
1H05C2	HBcAb Test Card	25 or 40
1H11C2	HBV-2 Panel Test Card (HBsAg/HBsAb)	25 or 40
1H10C2	HBV-3 Panel Test Card (HBsAg/HBeAg/HBcAb)	25
1H12C2	HBV-3 Panel Test Card (HBsAg/HBsAb/HBcAb)	25
1H0652	HBV-5 Panel Test Card	25
1H01C2-2	HBsAg Self Test Kit	1
1H0652-2	HBV-5 Self Test Kit	1

Dengue

One Step Dengue Test (Whole Blood / Serum / Plasma)

Rapid Dengue Test is a rapid immunochromatographic assay for the simultaneous detection of IgG and IgM antibodies to Dengue virus in human whole blood, serum or plasma. The product is used as a screening test for Dengue viral infection and an aid for differential diagnosis of primary and secondary infections in conjunction with other criteria.

Dengue IgM: Sensitivity: 95.8%, Specificity: 98.7%.

Dengue IgG: Sensitivity: 92.1%, Specificity: 98.6%.

Dengue NS1: Sensitivity: 96.1%, Specificity: 98.7%.

This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1N04S2	Dengue IgG/IgM Combo Test Strip	50 or 100
1N04C2	Dengue IgG/IgM Combo Test Card	25 or 40
1N23S2	Dengue NS1 Antigen Test Strip	50 or 100
1N23C2	Dengue NS1 Antigen Test Card	25 or 40
1N04C2-2	Dengue IgG/IgM Combo Self Test Kit	1
1N24C2-2	Dengue NS1 Antigen and IgG/IgM Duo Panel Self Test Kit	1

Rickettsia

One Step Rickettsia Test (Whole Blood / Serum / Plasma)

Rapid Rickettsia Test is a rapid immunochromatographic assay for the simultaneous detection of IgG and IgM antibodies to Rickettsia in human whole blood, serum or plasma. The assay is used as a screening test for Rickettsia infection and an aid for differential diagnosis of primary and secondary infections in conjunction with other criteria.

Sensitivity: 96%, Specificity: 99.3%.

This test is available in cassette format.



Cat. No.	Product Name	No. per Kit
1O06C2	Rickettsia IgG/IgM Combo Test Card	25 or 40

Malaria (pf, pv, pan)

One Step Malaria Antigen Test

Rapid Malaria Antigen Test is for the rapid qualitative determination of one or more of the known Malaria species: *P. falciparum*, *P. vivax*, *P. ovale*, and/or *P. malariae* by detecting lactate dehydrogenase (LDH) in human blood. This test is an aid in the diagnosis of Malaria infection.

Sensitivity: 95.1%, Specificity: 96.0%.

This test is available in cassette format.



Cat. No.	Product Name	No. per Kit
1Q05C2	Malaria pan Antigen Test Card	25 or 40
1Q06C2	Malaria pf Antigen Test Card	25 or 40
1Q07C2	Malaria pf/pv Antigen Test Card	25 or 40
1Q10C2	Malaria pf/pan Antigen Test Card	25 or 40
1Q06C2-2	Malaria pf Antigen Self Test Kit	1
1Q07C2-2	Malaria pf/pv Antigen Self Test Kit	1

H.Pylori (Antigen and Antibody)

One Step H.Pylori Antigen Test (Stool)

One Step H.Pylori Antibody Test (Whole Blood / Serum / Plasma)

Rapid H. pylori Antigen Test is an in vitro qualitative immunochromatographic assay for the rapid detection of *Helicobacter pylori* antigens in human stool specimen.


The test results are intended to aid in the diagnosis of H. pylori infection, to monitor the effectiveness of therapeutic treatment and confirm the eradication of H. pylori in peptic ulcer patients.

H.Pylori Antigen: Sensitivity: 94.0%, Specificity: 96.7%.

Rapid H. Pylori Antibody Test is a rapid immunochromatographic assay for the detection of antibodies to H. pylori in human serum, plasma or whole blood. The assay is to provide an aid in diagnosis of H. pylori infection.

H.Pylori Antibody: Sensitivity: 100%, Specificity: 97.0%.

This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1O02S1	H. Pylori Antibody Test Strip	50 or 100
1O02C1	H. Pylori Antibody Test Card	25 or 40
1O03C4	H.Pylori Antigen Test Card	25 or 40


Tuberculosis

One Step Tuberculosis Test (Whole Blood / Serum / Plasma)

Rapid Tuberculosis (TB) Test is a rapid chromatographic immunoassay for the qualitative detection of anti-TB (*M. tuberculosis*, *M. bovis* and *M. africanum*) antibodies (all isotypes: IgG, IgM, IgA, etc.) in human serum, plasma or whole blood.

Sensitivity: 97.4%, Specificity: 98.2%.

This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1O01S2	Tuberculosis Test Strip	50 or 100
1O01C2	Tuberculosis Test Card	25 or 40


Syphilis

One Step Syphilis Antibody Test (Whole Blood / Serum / Plasma)

Rapid Syphilis Antibody Test is a rapid immunochromatographic assay for the detection of antibodies to *Treponema pallidum* in human whole blood, serum or plasma. The assay is used as a screening test for *T. pallidum* infection (also known as Syphilis).

Sensitivity: 99.6%, Specificity: 99.93%.

This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1N02S2	Syphilis Antibody Test Strip	50 or 100
1N02C2	Syphilis Antibody Test Card	25 or 40
1N02C2-2	Syphilis Antibody Self Test Kit	1

Rotavirus

One Step Rotavirus Antigen Test (Stool)

Rapid Rotavirus Test is an in vitro qualitative immunochromatographic assay for the rapid detection of rotavirus antigens in human stool specimen. The test results are intended to aid in the diagnosis of rotavirus infection and to monitor the effectiveness of therapeutic treatment.

Sensitivity: 98%, Specificity: 97%.

This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1N07S4	Rotavirus Antigen Test Strip	50 or 100
1N07C4	Rotavirus Antigen Test Card	25 or 40


Adenovirus

One Step Adenovirus Antigen Test (Stool)

Rapid Adenovirus Test is an in vitro qualitative immunochromatographic assay for the rapid detection of adenovirus antigens in human stool specimen. The test results are intended to aid in the diagnosis of adenovirus infection and to monitor the effectiveness of therapeutic treatment.

Sensitivity: 96.6%, Specificity: 99.3%.

This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1N06C4	Adenovirus Antigen Test Card	25 or 40
1N08S4	Adeno-Rota Virus Antigen Combo Test Strip	50 or 100
1N08C4	Adeno-Rota Virus Antigen Combo Test Card	25 or 40

HAV

One Step Hepatitis A Antibody Test (Whole Blood / Serum / Plasma)

Rapid HAV IgG/IgM Combo Test is a single use, rapid device intended for qualitative detection of IgG-class and IgM-class antibodies to hepatitis A virus (HAV) in serum, plasma or whole blood samples. It is intended to be used in clinical laboratories for diagnosis of acute hepatitis A and management of patients related to infection with hepatitis A virus.

Sensitivity: 96.4%, Specificity: 99.3%.

This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1H15S2	HAV IgM Test Strip	50 or 100
1H15C2	HAV IgM Test Card	25 or 40
1H16S2	HAV IgG/IgM Combo Test Strip	50 or 100
1H16C2	HAV IgG/IgM Combo Test Card	25 or 40

HEV

One Step Hepatitis E Antibody Test (Whole Blood / Serum / Plasma)

Rapid HEV IgG/IgM Combo Test is a single use, rapid device intended for qualitative detection of IgG-class and IgM-class antibodies to hepatitis E virus (HEV) in serum, plasma or whole blood samples. It is intended to be used in clinical laboratories for diagnosis of acute hepatitis E and management of patients related to infection with hepatitis E virus.

Sensitivity: 100%, Specificity: 98.8%.


This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1H19S2	HEV IgM Test Strip	50 or 100
1H19C2	HEV IgM Test Card	25 or 40
1H20S2	HEV IgG/IgM Combo Test Strip	50 or 100
1H20C2	HEV IgG/IgM Combo Test Card	25 or 40

Other Infectious Diseases

Cat. No.	Product Name	No. per Kit
1N03C2	Hantavirus IgG/IgM Combo Test Card (HFRS)	25
1N17C5	Influenza A Test Card (Flu A)	20
1N18C5	Influenza A + Influenza B Test Card (Flu A+B)	20
1O04C2	S.Typhi IgG/IgM Duo Test Card	25
1O08C5	Strep A Test	25
1O09C2	S.Typhi Antigen Test Card	25
1O13C2	Salmonella typhi/paratyphi Antigen Duo Test	25
1O15S4	V.Cholerae O1 Test Strip	50 or 100
1O16S4	V.Cholerae O139 Test Strip	50 or 100
1Q08C1	Toxoplasma IgG/IgM Combo Test Card	25



Single (Cassette or Strip)



Cat. No.	Product Name	No. per Kit
1L01S3	Amphetamine Test Strip	50 or 100
1L01C3	Amphetamine Test Card	25 or 40
1L02S3	Barbiturate Test Strip	50 or 100
1L02C3	Barbiturate Test Card	25 or 40
1L03S3	Benzodiazepine Test Strip	50 or 100
1L03C3	Benzodiazepine Test Card	25 or 40
1L04S3	Buprenorphine Test Strip	50 or 100
1L04C3	Buprenorphine Test Card	25 or 40
1L05S3	Cocaine Test Strip	50 or 100
1L05C3	Cocaine Test Card	25 or 40
1L06S3	EDDP Test Strip	50 or 100
1L06C3	EDDP Test Card	25 or 40
1L07S3	Ketamine Test Strip	50 or 100
1L07C3	Ketamine Test Card	25 or 40
1L08S3	MDMA Test Strip	50 or 100
1L08C3	MDMA Test Card	25 or 40
1L09S3	Methadone Test Strip	50 or 100
1L09C3	Methadone Test Card	25 or 40
1L10S3	Methamphetamine Test Strip	50 or 100
1L10C3	Methamphetamine Test Card	25 or 40
1L11S3	Opiate Test Strip	50 or 100
1L11C3	Opiate Test Card	25 or 40
1L12S3	Opiate-2000 Test Strip	50 or 100
1L12C3	Opiate-2000 Test Card	25 or 40
1L13S3	Phencyclidine Test Strip	50 or 100
1L13C3	Phencyclidine Test Card	25 or 40
1L14S3	Marijuana (THC) Test Strip	50 or 100
1L14C3	Marijuana (THC) Test Card	25 or 40
1L15S3	TCA Test Strip	50 or 100
1L15C3	TCA Test Card	25 or 40
1L16S3	Oxycodone Test Strip	50 or 100
1L16C3	Oxycodone Test Card	25 or 40
1L30S3	Propoxyphene Test Strip	50 or 100
1L30C3	Propoxyphene Test Card	25 or 40
1L31S3	Tramadol Test Strip	50 or 100
1L31C3	Tramadol Test Card	25 or 40
1L35S3	Fentanyl Test Strip	50 or 100
1L35C3	Fentanyl Test Card	25 or 40
1L36S3	Methylphenidate Test Strip	50 or 100
1L36C3	Methylphenidate Test Card	25 or 40
1L44S3	Cotinine Test Strip	50 or 100
1L44C3	Cotinine Test Card	25 or 40
1L45S3	Clonazepam Test Strip	50 or 100
1L45C3	Clonazepam Test Card	25 or 40
1L46S3	Methaqualone Test Strip	50 or 100
1L46C3	Methaqualone Test Card	25 or 40
1L49S3	K2 Test Strip	50 or 100
1L49C3	K2 Test Card	25 or 40
1L29S5	Saliva Alcohol Test	25 or 100
1L29S3	Urine Alcohol Test Strip	50 or 100
1L29C3	Urine Alcohol Test Card	25 or 40

Saliva Test


Cat. No.	Product Name	No. per Kit
1L01C5	Amphetamine Saliva Test Card	10 or 25
1L03C5	Benzodiazepine Saliva Test Card	10 or 25
1L04C5	Buprenorphine Saliva Test Card	10 or 25
1L05C5	Cocaine Saliva Test Card	10 or 25
1L08C5	MDMA Saliva Test Card	10 or 25
1L09C5	Methadone Saliva Test Card	10 or 25
1L10C5	Methamphetamine Saliva Test Card	10 or 25
1L11C5	Opiate Saliva Test Card	10 or 25
1L14C5	Marijuana (THC) Saliva Test Card	10 or 25

More drugs of abuse saliva tests are available. Please contact us or your agents for details.



Multiple Drug Panel Test

Cat. No.	Product Name	No. per Kit
1L18M3	DOA-2 Panel Test, Dip Device	25 or 40
1L18C3	DOA-2 Panel Test, Cassette Device	25 or 40
1L19M3	DOA-3 Panel Test, Dip Device	25
1L19C3	DOA-3 Panel Test, Cassette Device	25
1L20M3	DOA-4 Panel Test, Dip Device	25
1L20C3	DOA-4 Panel Test, Cassette Device	25
1L21M3	DOA-5 Panel Test, Dip Device	25
1L21C3	DOA-5 Panel Test, Cassette Device	25
1L22M3	DOA-6 Panel Test, Dip Device	25
1L22C3	DOA-6 Panel Test, Cassette Device	25
1L23M3	DOA-7 Panel Test, Dip Device	25
1L23C3	DOA-7 Panel Test, Cassette Device	25
1L24M3	DOA-8 Panel Test, Dip Device	25
1L24C3	DOA-8 Panel Test, Cassette Device	25
1L25M3	DOA-9 Panel Test, Dip Device	25
1L25C3	DOA-9 Panel Test, Cassette Device	25
1L26M3	DOA-10 Panel Test, Dip Device	25
1L26C3	DOA-10 Panel Test, Cassette Device	25
1L33M3	DOA-11 Panel Test, Dip Device	25
1L34M3	DOA-12 Panel Test, Dip Device	25
1L28C3	Methamphetamine/Marijuana(THC) Test	25 or 40



DOA Test Cup

Drugs Of Abuse/Alcohol Test Cup, is an immunochromatography based one step in vitro test. It is designed for qualitative determination of drug substances in human urine specimens. This assay may be used in the point of care setting.

Cat. No.	Product Name	No. per Kit
1L27P3	DOA Panel Test Cup, up to 14 drugs, Split Type	1 cup/pouch
1L27K3	DOA Panel Test Cup, up to 14 drugs, On-site Type	1 cup/pouch



Cardiac

One Step Troponin I, Myoglobin, CK-MB Test (Whole Blood / Serum / Plasma)

Rapid Cardiac Panel Test is an immunochromatography based on one step in vitro test. It is designed for qualitative determination of cardiac troponin I (cTnI), CK-MB and Myoglobin in human serum, plasma or whole blood specimens as an aid in the diagnosis of myocardial infarction.

TnI:

Sensitivity: 98.4%, Specificity: 98%.

Myoglobin:

Sensitivity: 100%, Specificity: 99.3%.

CK-MB:

Sensitivity: 100%, Specificity: 99.5%.

This test is available in strip and cassette format.

Cat. No.	Product Name	No. per Kit
1J01C1	CK-MB Test Card - Serum	25 or 40
1J01S2	CK-MB Test Strip	50 or 100
1J01C2	CK-MB Test Card	25 or 40
1J02C1	Troponin I Test Card - Serum	25 or 40
1J02S2	Troponin I Test Strip	50 or 100
1J02C2	Troponin I Test Card	25 or 40
1J03C1	Myoglobin Test Card - Serum	25 or 40
1J03S2	Myoglobin Test Strip	50 or 100
1J03C2	Myoglobin Test Card	25 or 40
1J04C1	Cardiac Panel Test Card - Serum	25 or 40
1J04C2	Cardiac Panel Test Card	25 or 40



FOB

One Step Fecal Occult Blood Test (Fecal)

Rapid Fecal Occult Blood Test is a qualitative test that detects human hemoglobin in human fecal specimen. The test is a visual one step, in-vitro assay. It is intended for professional use to help diagnose gastrointestinal bleeding.

The analytical sensitivity of the test is 50 ng/mL h hemoglobin or 12.5 ug h hemoglobin/g feces

Sensitivity: 98.4%, Specificity: 100%.

This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1M01S4	Fecal Occult Blood Test Strip	50 or 100
1M01C4	Fecal Occult Blood Test Card,with 2-End Tube	25 or 40
1M01C4-2	Fecal Occult Blood Test Card,with 1-End Tube	25 or 40
1M01C4-3	Fecal Occult Blood Self Test Kit	1

PSA

One Step Prostate Specific Antigen Test (Whole Blood / Serum / Plasma)

Rapid PSA Test is an immunochromatography based on one step in vitro test. It is designed for the rapid semi-quantitative determination of human prostate specific antigen (PSA) in serum, plasma or whole blood specimens.

Sensitivity: 100%, Specificity: 98.3%.

This test is available in strip and cassette format.



Cat. No.	Product Name	No. per Kit
1C05S2	PSA Test Strip	50 or 100
1C05C2	PSA Test Card	25 or 40
1C05C2-2	PSA Self Test Kit	1


AFP

One Step Alpha-fetal Protein Test (Whole Blood / Serum / Plasma)

Rapid AFP Test is an immunochromatography based one step in vitro test. It is designed for the rapid qualitative determination of AFP in serum, plasma or whole blood specimens to aid in the diagnosis of Hepatocellular Carcinoma or fetal open neural tube defects.

Sensitivity: 99.4%, Specificity: 99.81%.

This test is available in cassette or Strip format.



Cat. No.	Product Name	No. per Kit
1C01S1	AFP Test Strip	50 or 100
1C01C2	AFP Test Card	25 or 40
1C01C2-2	AFP Self Test Kit	1


CEA

One Step Carcinoembryonic Antigen Test (Whole Blood / Serum / Plasma)

Rapid CEA Test is an immunochromatographic assay for the qualitative determination of Carcinoembryonic Antigen (CEA) in human serum, plasma or whole blood specimens. This device is intended as an aid in monitoring patients for disease progress or response to therapy or for the detection of recurrent or residual disease.

Sensitivity: 97.85%, Specificity: 99.52%.

This test is available in cassette or strip format.



Cat. No.	Product Name	No. per Kit
1C02S1	CEA Test Strip	50 or 100
1C02C2	CEA Test Card	25 or 40
1C02C2-2	CEA Self Test Kit	1

HCG

One Step HCG Pregnancy Test (Serum / Urine)

Rapid HCG Pregnancy Serum/Urine Test is a chromatographic immunoassay for the early detection of human chorionic gonadotropin (HCG) in human serum, plasma or urine specimens.

Sensitivity: 100%, Specificity: 100%.

This test is available in strip, cassette and midstream format.

Cat. No.	Product Name	No. per Kit
1D04S1	HCG Pregnancy Serum/Urine Test Strip	50 or 100
1D04C1	HCG Pregnancy Serum/Urine Test Card	25 or 40
1D04S3	HCG Pregnancy Urine Test Strip	50 or 100
1D04C3	HCG Pregnancy Urine Test Card	25 or 40
1D04M3	HCG Midstream Test	25 or 40
1D04C3-3	HCG Pregnancy Urine Self Test Kit-Cassette	1
1D04M3-3	HCG Pregnancy Urine Self Test Kit-Midstream	1

LH

One Step LH Ovulation Test (Urine)

Rapid LH Test is a rapid chromatographic immunoassay for the semi-quantitative detection of human luteinizing hormone (LH) in women's urine to predict the time of ovulation. The test is intended for healthcare professional and home use.

Sensitivity: 100%, Specificity: 100%.

This test is available in strip , cassette and midstream format.

Cat. No.	Product Name	No. per Kit
1D02S3	LH Test Strip	50 or 100
1D02C3	LH Test Card	25 or 40
1D02M3	LH Midstream Test	25 or 40
1D02C3-2	LH Self Test Kit-Cassette	1
1D02M3-2	LH Self Test Kit-Midstream	1



Questions? Please contact us at

info@bosonbio.com



Boson Biotech Co.,Ltd.

90-94 Tianfeng Road, Jimei North Industrial Park,
Xiamen, Fujian, China

Tel: 86-592-3965101

Fax: 86-592-3965155

E-mail: info@bosonbio.com

www.bosonbio.com